



Forever Medical

## YJ-U100

### Full Digital Palm Ultrasound Scanner



#### **Main Features:**

Full digital beam forming technology

Generate Report

Human and Veterinary double system(depend on your needs)

7 inch LCD LCD

A/DC and LI-battery operated

Light and compact design, simple to use

High quality keyboard and clear image

USB 2.0 for real-time picture uploading to PC

Probe automatic identification

Tel:+86-371-61315851

Email:[judy@yjcompany.com](mailto:judy@yjcompany.com)

Whatsapp/Skype/Phone:+86 15637181685



**Forever Medical**

Net weight: 1.2 Kgs

Dimension: 320×170×65mm

**Main Technical Specification:**

- Scanning Method: Convex, Micro-convex, Linear, Transvaginal
- Frequency for Convex : 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0
- Display Mode: B, B+B, B+M, M, 4B
- Gray Scale: 256 levels
- Gamma: 4
- Monitor: 7 inch LCD
- Operation interface: English/Chinese switchable
- Scanning Depth : 180mm
- Image magnification: \*1.0, \*1.2, \*1.3, \*1.5, \*1.6, \*1.8, \*2.0 times
- Image control: left/right, black/white, up/down
- General Measurement: Distance, Circumference/Area, Volume, Heart Rate
- 512 frames for cine-memory
- Human Body Mark : 15 kinds
- 4 focusing combination selectable
- Characteristic showing in real-time clock, ID, Sex, Age, Hospital name, Focusing Info and Gain (Near, Far, Overall) info.
- 8 Segments of TGC adjustment, full gain control
- Puncture function is available under B mode.
- Generate HistoG

**Human Measurement:**

OB : CRL

GS

BPD-FW

FL-FW

AC-FW

HC-FW



Forever Medical

APTD-FW

APD-FW

**Standard configuration:**

Main unit

Convex array probe

Li-battery

USB

**Options:**

Linear array probe

Micro-Convex array probe

Trans-vaginal array probe

Video Printer





Forever Medical



Tel:+86-371-61315851

Email:judy@yjcompany.com

Whatsapp/Skype/Phone:+86 15637181685



**Forever Medical**



**Tel:+86-371-61315851**

**Email:judy@yjcompany.com**

**Whatsapp/Skype/Phone:+86 15637181685**